

## Calibration of the CPS120

- Put Y<sub>2</sub>O<sub>3</sub> standard in the sample holder
  - Compress and level using a glass slide
  - Put sample holder into sample spinner. Make sure it slides all the way in
  - Close door
- To take the Pattern
  - In Acquisition
  - Clear the screen using the “clear” button
  - DO NOT pre-set the time (set time to 0sec) just watch the live time
  - Take pattern for 5mins: After 5mins (live time)
    - Open the door and remove the whole sample spinner (not just the sample holder)
    - Close door
    - Open Shutter
  - Take direct beam measurement for 45sec

Need to have a direct beam peak in same file as the pattern for calibration

- Save file as a DAT file (.dat)
  - Example:
    - 2007-16-05 caliby2o3xtalmirror.dat
      - Must save file with date (2007-16-05), type of file (calib), calibrant powder (y2o3) and what is in the path of the beam (xtal, mirror, or xtalmirror)
- Open the .dat file in Data Treatment (Imad)
- Find Peaks
  - In the “Fitting” tab
    - Peak Search
      - Use CF mode Put numbers 25 (far left box) and 6 (far right box) in the boxes to the left of the CF mode box
      - make sure fit channel and bck sub boxes are NOT checked
      - Set the range of 2theta to be from 0 to 120 in the boxes in the bottom right corner.
    - Click “Peak Search”
      - You need to make sure that all peaks are marked with a little upside down yellow triangle from 0-120°
        - To add peaks hold down the A (left click)
        - To remove peaks, hold down the Z (left click)
    - Under “View” click:
      - X-scale
      - channel

- Fit Peaks
  - Still in the “Fitting” tab
    - Choose Gaussian to fit
    - Make sure the boxes to the right of the Gaussian tab has the numbers 35 and 200 in them
    - Ensure that the following boxes are CHECKED
      - Fit Channel
      - Bck sub
  - Click on “Peak Fit”

Imad will do the peak fit and when done you will see dashed lines.

- Overlay of dashed lines (peak fitting positions) to solid lines (diffraction pattern)
  - Go to “VIEW” click
    - X-axis
    - Channel
  - Please note and remember the channel number of the direct beam, you will need it during calibration The value will be somewhere around the 320’s

- Calibration
  - Go to “Calibration” click
    - Cubic Spline

A new window will open: This is where you will create the calibration file

- In that new window select
  - YOBE/CU-D from the pull down menu on the bottom left
  - Set the channel value to the value of the direct beam that you noted earlier
  - Set “at angle (°)” to 0
  - Click the avg step carrots (<or>) to align the red arrows in the pattern with the blue arrows
    - You will never get all of the arrows to align up, so focus on aligning the arrows in the lower 2theta values (20-60°)
  - Watch the R-value (top right of window) You want/need a value below 0.1.
  - In the table, check the Standard and 2Tcalc columns to make sure the values are similar
  - Confirm the wavelength is correct in the bottom left hand corner of the screen. This value will be used in the subsequent generation of the calibration file:
    - Cu  $K_{\alpha 1}$  (if with monochromator) = 1.540598 Å
    - Cu  $K_{\alpha, avg}$  (without monochromator) = 1.54184 Å

- When you are satisfied
  - Click on the “Create Cal” button

This creates a .cal file with the same name as the .dat file

- Double-Checking
  - Close the files that you are working with in Imad
    - File
    - Close
  - In Imad Open the file that you just closed

You should see red dashes on it now

- In the “standard” tab
  - Under the “Standard” heading
    - Select the calibrant powder (YOBE/CU-D)

It will overlay, and you will see blue dashes that will not line up with your peaks

- Go into “Edit”
  - Change Cal
  - Select the .cal file that you just made
    - Remember it will have the same name as the .dat file

Now the peaks in the diffraction pattern will line up with the blue dashed lines

- To Save
  - “Save File”
  - Chose DAT file
  - Overwrite the .dat file that you made earlier

You do this so the .dat file is calibrated.

- In Acquisition
  - Click on the “Change Calibration” button in the top left corner
  - Select the .cal file that was just made
  - Close Acquisition
  - Open Acquisition
    - Double check that the calibration file that is displayed is correct
- Put the sample spinner back in place
- **BOLDLY** write the name of the calibration file on the record sheet